









SKIN PRODUCTS

PRODUCT CODE	PRODUCT DESCRIPTION	PRODUCT
LITA18/H050	Diphoterine® Solution MICRO 100ml Spray <ul style="list-style-type: none"> Recommended for area size of a hand or face. Begin rinsing within 60 seconds of chemical splash. Use entire contents of package. 	
LITA18/H051	Diphoterine® Solution MINI 200ml Spray <ul style="list-style-type: none"> Recommended for area size of a limb or front torso. Begin rinsing within 60 seconds of chemical splash. Use entire contents of package. 	
LITA18/H057	Diphoterine® Solution DAPD 5Lt Portable Shower <ul style="list-style-type: none"> Recommended for full body decontamination. Begin rinsing within 60 seconds of chemical splash. Use entire contents of package. 	
LITA18/H069	Diphoterine® Exchange Serviced Portable Shower	


EYE PRODUCTS

PRODUCT CODE	PRODUCT DESCRIPTION	PRODUCT
LITA18/H055	Diphoterine® Solution LIS 50ml Individual Eyewash <ul style="list-style-type: none"> Begin rinsing within 10 seconds of chemical splash. Use entire contents of package. 	
LITA18/H058	Diphoterine® Solution LPD 500ml Portable Bag with Eye Cup <ul style="list-style-type: none"> Begin rinsing within 60 seconds of chemical splash. Use entire contents of package. 	
LITA18/H058	Diphoterine® Solution LPM 500ml Eyewash Bottle <ul style="list-style-type: none"> Begin rinsing within 60 seconds of chemical splash. Use entire contents of package. 	

EYE & SKIN WASH STATION

PRODUCT CODE	PRODUCT DESCRIPTION	PRODUCT
LITA18/H103	Diphoterine® Solution LPMD. DAP Wall Mount Wash Station <ul style="list-style-type: none"> 2 x 500ml Eyewash Bottles. 1 x 200ml Mini Spray. 	
	Diphoterine® All in One Case DAP. LPMD. DAP <ul style="list-style-type: none"> A portable shower (DAP 5 litres) for full body rinsing. Two 500 ml bottles (LPMD). One 200 ml MINI DAP. One bottle of Afterwash II® solution. 	

POUCHES

PRODUCT CODE	PRODUCT DESCRIPTION	PRODUCT
LITA 18/H075	Pouch Hi viz to suit 200ml Diphoterine Spray	
LITA 18/H073	Diphoterine® Belt 38mm to suit	